



NexxZr® Shine

高效抛光。
自然光泽。

高效抛光。
自然光泽。

NexxZr® Shine



NexxZr. Shine是一种两步法抛光系统用于有效抛光牙科陶瓷,例如:氧化锆、氧化铝、二硅酸锂或玻璃陶瓷。

混入了相同尺寸的金刚石的橡胶抛光工具,工具之间配合高效,从而得到了用户的信任。此外它们具有较长的使用寿命,并且可以为您的修复体展现栩栩如生的自然光泽。

NexxZr Shine 使用在整体和部分贴面修复体均匀的表面。三种优化后不同结构和形状抛光工具分别为: Turn 18、Twist 17 和 Tip 5。抛光工具提供光滑(第1步,预抛光)和光泽(第2步,高光抛光)的颗粒度。通用的金刚石抛光膏 NexxZr Shine Paste 完善了产品系列。它可以去除微划痕并留下光滑如玻璃的表面,具有高光泽。

技术数据



类型
两步金刚石抛
光系统



材料
橡胶内混入金刚石
抛光工具



颗粒细度
柔滑 (smooth),
光滑 (gloss)



适用于
所有陶瓷



3个形状,1个抛光膏
Turn 18, Twist 17,
Tip 5, 抛光膏

1 两步金刚石抛光系统 轻松高效地可用于所有陶瓷材料

步骤 1
预抛光



颗粒度：
柔滑 (smooth)



外观：
丝光



步骤 2
高光抛光



颗粒度：
光滑 (gloss)



外观：
高光泽



2 橡胶工具内混入尺寸相等的金刚石 栩栩如生的光泽,使用寿命长



镜面般光滑、自然的光泽效果

3 三种不同形状的抛光工具和一个抛光膏 抛光效果佳,不受结构和表面限制



Turn圆盘18
用于曲面



Twist 刷轮 17
用于表面和沟窝



Tip圆锥5
用于表面和沟窝



抛光膏
用于表面和沟窝



交货清单



套装

NexxZr Shine Starter Kit
Turn 18 smooth / 1pcs.
Turn 18 gloss / 1pcs.
Twist 17 smooth / 2pcs.
Twist 17 gloss / 2pcs.
Tip 5 smooth / 2pcs.
Tip 5 gloss / 2pcs.
Paste 20 g / 1pcs.

NexxZr Shine Polishing Kit
Turn 18 smooth / 1pcs.
Turn 18 gloss / 1pcs.
Twist 17 smooth / 2pcs.
Twist 17 gloss / 2pcs.
Tip 5 smooth / 2pcs.
Tip 5 gloss / 2pcs.



补充装

NexxZr Shine Turn 18 smooth / 5pcs.
NexxZr Shine Turn 18 gloss / 5pcs.
NexxZr Shine Twist 17 smooth / 5pcs.
NexxZr Shine Twist 17 gloss / 5pcs.
NexxZr Shine Tip 5 smooth / 5pcs.
NexxZr Shine Tip 5 gloss / 5pcs.
NexxZr Shine Paste 20 g / 1pcs.

Manufacturer

Sagemax Bioceramics, Inc.
34210 9th Ave. South, Suite 118
Federal Way, WA 98003, USA
sagemax.com